HyZone

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

HyZone Class II Microbiological Safety Cabinets

HyZone Class II Microbiological Safety Cabinets from ThermoLab blend knowledge, experience and innovative design changes in creating state-of-the-art equipment as part of its offerings. HyZone is a Class II Microbiological Cabinet, designed and built to meet, if not exceed, performance requirements of the EN 12469:2000 Standard. 70% air is re-circulated through the main H14 HEPA Filter within the cabinet, while the remaining 30% is exhausted through two H14 HEPA filters in series.

HyZone Cabinets are suitable for handling microorganisms and pathogens as defined by the appropriate International Standards, current health and safety guidelines and legislation aimed at safeguarding health and safety of operators at work.

95 % of your requirements are met by Class II A2 Cabinets. However, Class I and Class III cabinets can be manufactured if your work so warrants.

Salient features of HyZone Class II Microbiological Safety Cabinets


- Independent blowers for Supply and Exhaust air
- Double HEPA filtration of Exhaust air
- Negative pressure envelope (shroud) for HEPA filter Plenums and Return air path to prevent escape of Contaminated air to atmosphere and work area
- Simultaneous real time monitoring and display of Supply and Exhaust air velocities
- Monitoring of work area temperature
- Sliding front sash with counter weights for easy operation
- Microprocessor controlled operations with Normal, Fast Start and Service modes and alarms
- Ergonomic design like angled front sash, arm rest etc., for operator comfort
- RS 232 connectivity
- EN 12469:2000 compliant

Soft Touch Keyboard

Membrane touch keyboard and the rear-lit LCD display all required data keeping the operator constantly up to date of the cabinet conditions:

- Estimated life of HEPA filters
- Estimated life of UV lamp
- Cabinet Temperature
- Laminar airflow (Unidirectional Down Flow) velocity and frontal air barrier (in flow) velocity
- Front sash position

Audio Visual alarms are provided for :

- Out of range or incorrect laminar airflow velocity and frontal air barrier velocity
- Incorrect position of front sash-window
- Clogging of HEPA filters
- End of life-cycle of UV Lamp
- Out of range temperature

APPLICATIONS AREAS :

- Microbiology
- Medical Devices
- Virology
- Pharmaceuticals
- Haematology
- Pathology
- Cell culture
- Blood Bank
- Recombinant – DNA
- Eye Bank
- Stem Cells Research
- Reconstitution of Parenteral drugs
- Biotechnology
- Genetics.

Quality Assurance

Each unit is factory tested to meet international guidelines and dispatched in ready-to-plug-in condition. ThermoLab's testing services are available to assist you in routine maintenance and revalidation.

Corporate Office: Thermolab House, Plot No. 19, Vasai Municipal Ind. Area, Umela Road,
Vasai (West), Dist.-Thane, Maharashtra, India.

Registered Office: 1-A, Ojha Niwas, Patel Estate Marg, Jogeshwari (W), Mumbai-400102.
(E) info@thermolabscientific.com | (T) +91-250-2323156, 2324866/7 | (F) +91-250-2321656
www.thermolabscientific.com | www.thermolabgroup.com

"As part of our continuous improvement program, design specifications and photographs given in this brochure are subject to changes without notice and without attracting obligations for free modifications to equipment previously sold"